CoHere

Webinar

Boosting WASPLab® with
PhenoMATRIX®

Microbiology labs experience with A.I.

11 November 2021
16.00 to 17.00 CET


**PhenoMATRIX:
real-life
experience
of a high volume
bacteriology lab**

Simone Ambretti

Unit of Clinical Microbiology

**S.Orsola
University Hospital,
Bologna, Italy**

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
PHENOMATRIX : REAL LIFE EXPERIENCE OF A HIGH-VOLUME BACTERIOLOGY LAB

Unit of Clinical Microbiology S.Orsola University Hospital, Bologna, Italy

- ✓ **Microbiology lab for all Bologna metropolitan area (northern Italy, 1 million inhabitants)**
- ✓ **13 different hospitals**
- ✓ **Total of almost 4.000 beds**
- ✓ **Open 7/7, from 7.30 a.m. to 8 p.m.**
- ✓ **2020 diagnostic activity**
 - **140.000 blood cultures**
 - **100.000 CRE surveillance cultures**
 - **75.000 urine cultures**
 - **20.000 yeast cultures**
 - **7.000 GBS surveillance cultures**


Unit of Clinical Microbiology, S.Orsola University Hospital, Bologna, Italy

OUR JOURNEY TO AUTOMATION IN BACTERIOLOGY

1° step: AUTOMATION OF STREAKING

➤ 2009 **WASP** (Walk-Away Specimen Processor)

from manual to automated streaking of liquid based samples (urine, stool, swabs) + enrichment broth inoculation + slide preparation

Liquid Based Microbiology

WASP

OUR JOURNEY TO AUTOMATION IN BACTERIOLOGY

2° step: AUTOMATION OF INCUBATION AND PLATE READING

➤ **2016 WASPLAB**

OUR JOURNEY TO AUTOMATION IN BACTERIOLOGY

3° step: AUTOMATION OF PLATE INTERPRETATION

Digital
microbiology

- 2017 Implementation of **segregation protocols** (NEGATIVE vs POSITIVE) for CRE surveillance rectal swabs and urine cultures

Contents lists available at ScienceDirect

Journal of Global Antimicrobial Resistance

journal homepage: www.elsevier.com/locate/jgar

Rectal screening for carbapenemase-producing Enterobacteriaceae: a proposed workflow

Claudio Foschi^{a,b,*}, Paolo Gaibani^b, Donatella Lombardo^b, Maria Carla Re^{a,b},
Simone Ambretti^b

^a Microbiology Unit, DIMES, University of Bologna, via Massarenti 9, Bologna, Italy

^b Microbiology Unit, S. Orsola-Malpighi Hospital, Via Massarenti 9, Bologna, Italy

ARTICLE INFO

Article history:

Received 19 August 2019

Received in revised form 10 October 2019

Accepted 14 October 2019

Available online 19 October 2019

Keywords:

CPE

Carbapenemase

Klebsiella pneumoniae

MALDI-TOF/MS

WASPLab

Rectal screening

ABSTRACT

Objectives: Active screening is a crucial element for the prevention of carbapenemase-producing Enterobacteriaceae (CPE) transmission in healthcare settings. Here we propose a culture-based protocol for rectal swab CPE screening that combines CPE detection with identification of the carbapenemase type.

Methods: The workflow integrates an automatic digital analysis of selective chromogenic media (WASPLab[®]; Copan), with subsequent rapid tests for the confirmation of carbapenemase production [i.e. detection of *Klebsiella pneumoniae* carbapenemase (KPC)-specific peak by matrix-assisted laser desorption/ionisation time-of-flight mass spectrometry (MALDI-TOF/MS) or a multiplex immunochromatographic assay identifying the five commonest carbapenemase types]. To evaluate the performance of this protocol in depth, data for 21 162 rectal swabs submitted for CPE screening to the Microbiology Unit of S. Orsola-Malpighi Hospital (Bologna, Italy) were analysed.

Results: Considering its ability to correctly segregate plates with/without Enterobacteriaceae, WASPLab Image Analysis Software showed globally a sensitivity and specificity of 100% and 79.4%, respectively. Of the plates with bacterial growth ($n = 901$), 693 (76.9%) were found to be positive for CPE by MALDI-TOF/MS (KPC-specific peak for *K. pneumoniae*) or by immunochromatographic assay. Only 2.8% (16/570) of KPC-positive *K. pneumoniae* strains were missed by the specific MALDI-TOF/MS algorithm, being detected by the immunochromatographic assay. The mean turnaround time needed from sample arrival to the final report ranged between 18 and 24 h, representing a significant time saving compared with manual reading.

Conclusion: This workflow proved to be fast and reliable, being particularly suitable for areas endemic for KPC-producing *K. pneumoniae* and for high-throughput laboratories.

© 2019 International Society for Antimicrobial Chemotherapy. Published by Elsevier Ltd. This is an open access article under the CC BY-NC-ND license (<http://creativecommons.org/licenses/by-nc-nd/4.0/>).

PhenoMATRIX: real-life experience of a high volume bacteriology lab


Fig. 1. Workflow for rectal carbapenemase-producing Enterobacteriaceae (CPE) screening.

PHENOMATRIX : REAL LIFE EXPERIENCE OF A HIGH-VOLUME BACTERIOLOGY LAB


- 2020 introduction of **PhenoMatrix**: development of different protocols of automated plate-reading (CRE surveillance, GBS surveillance, yeast cultures, urine cultures)

**ARTIFICIAL
INTELLIGENCE
INTERPRETATION**

Algorithms

Lis Integration


Applied Artificial Intelligence

AUTOMATIC INTERPRETATION (PRESUMPTIVE RESULT)


PHENOMATRIX : REAL LIFE EXPERIENCE OF A HIGH-VOLUME BACTERIOLOGY LAB

PHENOMATRIX PROTOCOLS




✓ CRE surveillance rectal swabs

✓ GBS surveillance vagino-rectal swabs

✓ Yeast cultures
(genital and respiratory samples)

✓ Urine cultures


CRE surveillance rectal swabs

2020: 100.000 samples
Daily mean: > 330 swabs / day
Peaks: >800 swabs /day
Positivity rate: <5%

CULTURE PROTOCOL

Material = Rectal swab
Analysis = CPE screening

VACUTEST KIMA®, Chromagar KPC

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
ANSA: 10 µl

INCUBATORE: O₂

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IMAGE ANALYSIS ALGORITHM


- ✓ Qualitative evaluation
- ✓ Colony detection based on size and colour

- Pink/red (*Escherichia coli*)
- Blue/green (KESC)
- Transparent/white (suspected for *Acinetobacter baumannii*) **BUT ONLY FOR ICU PATIENTS**

PHENOMATRIX : REAL LIFE EXPERIENCE OF A HIGH-VOLUME BACTERIOLOGY LAB

Samples segregated into 3 different folders: **NEGATIVE**, **CHROM-ID**, **MIXED**

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
Negativa


- Batch validation of plates results
- Negative results sent to LIS
- Plates sent to waste

PHENOMATRIX : REAL LIFE EXPERIENCE OF A HIGH-VOLUME BACTERIOLOGY LAB

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- Batch validation of plates results
- **positive result and presumptive ID sent to LIS**
- plates sorted to follow-up activities column

PHENOMATRIX : REAL LIFE EXPERIENCE OF A HIGH-VOLUME BACTERIOLOGY LAB

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PHENOMATRIX : REAL LIFE EXPERIENCE OF A HIGH-VOLUME BACTERIOLOGY LAB

GBS surveillance recto-vaginal swabs


2020: 7.000 samples
Daily mean: 25 swabs / day
Positivity rate: 16-18%

CULTURE PROTOCOL

Material = Vagino-Rectal swab

Analysis = GBS screening

CHROMagar™ StrepB, VACUTEST KIMA

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ANSA: 1 microlitro


INCUBATORE: ossigeno

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IMAGE ANALYSIS ALGORITHM

- ✓ Qualitative evaluation
- ✓ Colony detection based on size and colour

Purple colonies =
suspected for *Streptococcus agalactiae*)


Performance of PhenoMatrix for the detection of Group B *Streptococcus* from recto-vaginal swabs

Claudio Foschi^{a,b,*}, Gabriele Turello^a, Tiziana Lazzarotto^{a,b}, Simone Ambretti^b

^a Microbiology, DIMES, University of Bologna, Bologna, Italy

^b Microbiology Unit, IRCCS S.Orsola-Malpighi Hospital, Bologna, Italy

Digital images of plates were automatically screened by Pheno-Matrix.

This software analyses the plates to identify differences in growth and colony colour.

A specific protocol was developed evaluating the medium type and the incubation time used by the laboratory software automatically sorted the cultures into the categories 'negative for GBS' and 'potential-positive for GBS'.

In our protocol, plates with no bacterial growth or with blue/colourless colonies are marked as GBS-negative, whereas plates with pink/light red/purple colonies are defined as presumptive GBS-positive.

PHENOMATRIX : REAL LIFE EXPERIENCE OF A HIGH-VOLUME BACTERIOLOGY LAB

2

C. Foschi et al. / Diagnostic Microbiology and Infectious Disease 101 (2021) 115427


Fig. 1. . Images of representative cultures on CHROMagar StrepB medium. (A-B) Streptococcus agalactiae - GBS; (C) Streptococcus anginosus; (D) Enterococcus faecalis; (E) Streptococcus oralis; (F) Staphylococcus epidermidis.

Cultures were also reviewed manually by a technologist to check software analyses and scored for the presence or absence of colonies resembling GBS.

Colonies consistent with GBS were confirmed by matrix-assisted laser desorption ionization–time of flight (MALDI-TOF) identification.

A specimen was considered a true-positive if a colony was confirmed to be GBS by MALDI-TOF.

A total of 587 plates (55%) were correctly segregated as GBS-negative by PhenoMatrix, with **no false-negative results**. The automatic reading system marked 481 plates (45%) as potentially positive; 158 of them were confirmed as GBS by MALDI-TOF.

The remaining cases (323; 30.2% of the total) sorted positive by PhenoMatrix were characterized by the absence of pink colonies at visual inspection (212 cases) or by a MALDI-TOF identification different from GBS (111 cases).

YEAST CULTURES: swabs, genital and respiratory samples

2020: 20.000 samples
Daily mean: 65 swabs / day


CULTURE PROTOCOL

Material = swabs, genital and
respiratory samples

Analysis = Yeast culture

IMAGE ANALYSIS ALGORITHM

✓ Quantitative evaluation


Colony detection based on colour:

- Green = *Candida albicans*
- blue = *Candida tropicalis*
- white/pink/purple = Other *Candida* species

Colonies counting:

- Different cut-off values for different types of samples

VACUTEST KIMA®, Chromagar Candida

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ANSA: 10 µl

INCUBATORE: O₂

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PHENOMATRIX : REAL LIFE EXPERIENCE OF A HIGH-VOLUME BACTERIOLOGY LAB

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-NO GROWTH OR
-*Candida* colonies count lower than
sample type- specific cut-off

- Batch validation of plates results
- negative result sent to LIS
- plates sent to waste

Introduction of Phenomatrix protocols for plate-reading automation of CRE and GBS Screening Specimens, Urine Cultures and Yeast Cultures on chromogenic media: real life experience of an high volume bacteriology lab

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GROWTH OF *Candida albicans* or *Candida tropicalis*
with colonies count higher than sample type- specific cut-off

- Batch validation of plates results
- positive result, colonies count, and species ID sent to LIS
- plates sorted to back-up activities column

GROWTH OF *Candida sp.* different from *Candida albicans* or *Candida tropicalis*
with colonies count higher than sample type- specific cut-off

- Batch validation of plates results
- positive result, colonies count, and presumptive ID (*Candida sp.*) sent to LIS
- plates sorted to MALDI-TOF ID column


URINE CULTURES

2020: 75.000 samples
Daily mean: 250 swabs / day

CULTURE PROTOCOL

Material = urine

Analysis = urine culture

VACUTEST KIMA®, Chromagar Orientation

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
INCUBATORE: O₂

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IMAGE ANALYSIS ALGORITHM

- ✓ **Quantitative evaluation**
- ✓ **High level of complexity**

- o *Escherichia coli*
- o *Enterococcus*
- o KESC
- o PMP
- o *Staphylococcus saprophyticus*
- o *Streptococcus agalactiae*
- o *Pseudomonas*
- o colonie bianche
- o colonie trasparenti
- o microcolonie
- o sconosciuto

Result evaluation based on:

- **Global colonies counting**
- **Number of different types of isolates**
- **Count, colour and size evaluation of single type of colonies**
- **Sex and age of patients (for GBS)**

NOTE: LIS limitation did not allow to use other types of parameters (leucocytes count)

URINE CULTURES


RESULT EVALUATION BASED ON COLONIES PREVALENCE RULES

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URINE CULTURES

RESULT EVALUATION BASED ON PATIENTS DATA

EVALUATE AS POSITIVE FOR GBS

With lower colonies count (10^4 CFU/ml)

ONLY IF

Sex: Female

Age: >14, <50 YEARS

PHENOMATRIX : REAL LIFE EXPERIENCE OF A HIGH-VOLUME BACTERIOLOGY LAB

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GROWTH OF MONOMICROBIAL
Escherichia coli or other gram-negatives
with colonies count higher than sample type- specific cut-off

- Batch validation of plates results
- positive result, colonies count, and species ID or presumptive ID sent to LIS
- plates sorted to follow-up activities column (AST or MALDI-TOF ID plus AST)


GROWTH OF MONOMICROBIAL
Enterococci or Staphylococci (*saprophyticus* or *aureus*)
with colonies count higher than sample type- specific cut-off

- Batch validation of plates results
- positive result, colonies count, and presumptive ID sent to LIS
- plates sorted to follow-up activities column (MALDI-TOF ID plus AST)

PHENOMATRIX : REAL LIFE EXPERIENCE OF A HIGH-VOLUME BACTERIOLOGY LAB

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
**POLYMICROBIAL GROWTH WITH THREE OR MORE
DIFFERENT TYPE OF COLONIES**

- Batch validation of plates results
- **CONTAMINATION** result sent to LIS
- plates sorted to back-up activities column

PHENOMATRIX : REAL LIFE EXPERIENCE OF A HIGH-VOLUME BACTERIOLOGY LAB

URINE CULTURES

LETTURA

Regole generali:

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
1 isolato,
stafilo/entero/coli/bacneg
10-39

2 isolati > 39 (totale)

2 isolati (no Strepto B) 10-39

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**MONOMICROBIAL GROWTH OF POSSIBLE
PATHOGENS WITH INTERMEDIATE TOTAL COUNT**

OR

**POLYMICROBIAL GROWTH WITH TWO DIFFERENT
TYPES OF COLONIES**

- Batch validation of plates results

- Plates sent to operator visual reading

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


Visualizza 10 elementi

More 80% of WASPLAB samples are currently evaluated using a Phenomatrix plate-reading automated protocol

OUR JOURNEY TO AUTOMATION IN BACTERIOLOGY

➤ Urine:
WASP+WASPLAB+PHENOMATRIX

➤ Rectal swabs:
WASP+WASPLAB+PHENOMATRIX

➤ Yeast and SGB cultures:
WASP+WASPLAB+PHENOMATRIX

➤ Respiratory samples: WASP (streaking +
slide preparation) + WASPLAB


➤ Other swabs and genital samples:
WASP+WASPLAB

➤ Positive blood cultures: **WASP** (streaking
+ slide preparation)


➤ Stool cultures: **WASP** (streaking and
broth inoculation)

➤ Other samples: manual streaking and
visual reading (traditional)

Sample distribution (2020 data)


Level of automation

WHERE WE WERE, WHERE WE ARE, WHERE WE ARE GOING?

- Over the last decade we have implemented **FULL LAB AUTOMATION** in our bacteriology lab
- Starting from LBM, we introduced automated streaking (**WASP**), automated incubation and plate reading (**WASPLAB**) and automated sample interpretation (**PHENOMATRIX**)
- Currently almost 80% of bacteriology samples are routinely evaluated by AI protocols
- Is our journey into automation over? No!!!! Next steps:
 - Improvement of existing protocols
 - Development of new protocols (multiple plates samples)

AUTOMATION = STANDARDIZATION, TRACEABILITY, EFFICIENCY

CoHere

Webinar

Boosting WASPLab® with
PhenoMATRIX®

Microbiology labs experience with A.I.

**THANK
YOU FOR
YOUR
ATTENTION**

BOLOGNA